4TEK SRL
Viale Dell'industria 66
20037 Paderno Dugnano (MI) ITALY
TEL +39 02 91082211
Fax +39 02 91084668
stefano@4-tek.it gabriella@4-tek.it riccardo@4-tek.it
www.4-tek.it

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UNITS (with cabinet)
UNITA' (caremate)
A004/L - A004
A005
A006/L - A006
OSM/L - OSM

<p>Contents</p> <p>SECTION A GENERAL FEATURES</p> <p>1 Description</p> <p>2 Type of appliances and selection</p> <p>3 System's features</p> <p>3.1 Principles of operation</p> <p>3.2 Installation criteria</p> <p>4 Installation design</p> <p>4.1 Aspiration circuit</p> <p>4.2 Liquid waste</p> <p>4.3 Air exhaust</p> <p>4.4 Electrical wiring</p> <p>SECTION B INSTALLATION USE AND MAINTENANCE</p> <p>5 General</p> <p>6 Requirements</p> <p>6.1 Use</p> <p>6.2 Aspiration pressure</p> <p>6.3 Ambient temperature</p> <p>7 Installation</p> <p>7.1 Transportation</p> <p>7.2 Assembly</p> <p>7.3 Electrical connection</p> <p>7.4 Starting</p> <p>8 Maintenance</p> <p>8.1 Maintenance by the user</p> <p>8.2 Maintenance by qualified personnel</p> <p>9 Troubleshooting</p> <p>10 Guarantee</p> <p>SECTION C SPECIAL INSTRUCTIONS ONLY FOR A006L</p> <p>11 General features</p> <p>12 Installation</p> <p>13 Operation</p> <p>14 Maintenance</p> <p>15 Control board</p> <p>_ Table 2 – Troubleshooting</p> <p>SECTION D SAFETY INSTRUCTIONS</p> <p>16 Introduction</p> <p>16.1 Foreword</p> <p>16.2 Definitions</p> <p>16.3 Symbols</p> <p>16.4 Qualified personnel</p>	<p>17 General</p> <p>17.1 Safety notes</p> <p>17.2 Use</p> <p>18 Installation</p> <p>18.1 General safety rules</p> <p>18.2 Installation</p> <p>18.3 Electrical connection</p> <p>18.4 Water and waste lines</p> <p>18.5 Starting</p> <p>19 Operation</p> <p>19.1 General safety rules</p> <p>20 Maintenance and repairs</p> <p>20.1 User's maintenance</p> <p>20.2 General safety rules</p> <p>20.3 Dismanting</p> <p>20.4 Reassembling</p> <p>- Table 3 – Technical Data</p> <p>- Table 4 – Symbols used on the equipment</p> <hr/> <p>SECTION A GENERAL FEATURES</p> <p>1 DESCRIPTION</p> <p>The Air Ring aspiration system comprises appliances, Aspirators and Units, the main component of which is a lateral channel blower with a radial blades impeller.</p> <p>The vacuum generated by the blower produces a depression that sucks secretions via the cannula in the dental unit.</p> <p>This economic system features simple construction, easy operation and reduced maintenance.</p> <p>2 TYPE OF APPLIANCE AND SELECTION</p> <p>The appliances covered by this manual are: Aspirators with and without cabinet.</p> <p>Table 1 gives the number of surgeries that can be operated simultaneously depending on the Model of appliance Installed</p> <table border="1"> <thead> <tr> <th rowspan="2">Model</th> <th rowspan="2">Surgeries simultaneous operation</th> <th colspan="2">Pipes</th> </tr> <tr> <th>Main P</th> <th>Secondary D</th> </tr> </thead> <tbody> <tr> <td>A001/L or A001 A004/ or A004 A006L or A006 OSM/L or OSM</td> <td>1</td> <td>mm 32</td> <td>mm 32</td> </tr> <tr> <td>A002 A005</td> <td>2</td> <td>50</td> <td>40</td> </tr> <tr> <td>A003</td> <td>3</td> <td>60</td> <td>50</td> </tr> <tr> <td>A007</td> <td>4-6</td> <td>80</td> <td>72</td> </tr> </tbody> </table> <p>Table 1</p>	Model	Surgeries simultaneous operation	Pipes		Main P	Secondary D	A001/L or A001 A004/ or A004 A006L or A006 OSM/L or OSM	1	mm 32	mm 32	A002 A005	2	50	40	A003	3	60	50	A007	4-6	80	72
Model	Surgeries simultaneous operation			Pipes																			
		Main P	Secondary D																				
A001/L or A001 A004/ or A004 A006L or A006 OSM/L or OSM	1	mm 32	mm 32																				
A002 A005	2	50	40																				
A003	3	60	50																				
A007	4-6	80	72																				

3 SYSTEM'S FEATURES

3.1 Principles of operation

Pic.1 shows the main components of the system and their respective function. Once the general switch on the mains has been set to ON, by simply lifting or replacing a tip in its stand, the operator actuates the low tension micro-switch (situated in the tip-holder) that starts or stops the aspiration. Such an arrangement allows the pump to run only for the required time with a consequent saving in power. The pump designed to stand heavy duties and for continuous or discontinuous running.

3.2 Installation criteria

1. Aspirators with cabinet can be installed in the surgery itself, when a better location is not available and are normally placed in the cabinet under the wash stand. In this case it is important to make sure that a proper ventilation is guaranteed. The pump shall not be bolted onto the floor. Its antivibration rubber pads ensure a smooth and very low noise operation. For those instances when an absolutely noiseless operation has to be achieved, an additional silencer (Art. S203) can be supplied to special order.

2. Having regards to their noise level, aspirators without cabinet are to be installed in a location closest to the surgery.

3. The area or location where the appliances are to be installed should feature:

- power mains adequate for the load indicated on the pump or pump plate
- A low tension control line from the surgery
- An adequate arrangement for the air exhaust
- An adequate waste line to deal with secretions, should also be provided (see Par 4.2)

4. Warning! In order to ensure proper ventilation of the blower and to avoid any risk of overheating or burning it is mandatory that an over-pressure valve be installed for machines having a power rating over 0,75 kW.

4 INSTALLATION DESIGN

4.1 Aspiration circuit

The diameter of suction pipes should be consistent with the air flow value. The layout of the waste lines, if required according to par 4.2 below, should be such as not to hamper the normal flow of liquids.

The main pipe P (Tab. 1) acts as a manifold and has different diameters in proportion to the volume of fluids coming from the dental units. The secondary pipes D that connect the dental units to the main pipe have smaller diameters.

It is advisable to install the piping with a slight slope towards the pump. Suitable inspection traps may be provided at selected locations.

It is also advisable to install the appliance in a lower position with respect to the suction pipes in order to improve drainage and consequently, its operation.


6.4 Liquid waste

(see Pic.1)

All air ring aspiration appliances require that an air/liquids separator **VS** must be installed between the dental unit **C** and the aspirator **A**.

In the case of A006L or A006 a separator is a component of the aspirator unit.

In any case the integration of this separator in the system is mandatory in order that no liquid would reach the aspiration circuit. Consequently it is necessary that a waste line **L** be foreseen near the separator or the Unit.

6.5 Air exhaust

In Units with cabinets air is exhausted through a silencer which is usually assembled under the blower.

6.6 Electrical wiring

The single or three-phase power line (as the case may be) must be consistent with the power of the Unit(s) installed. All national or local regulations concerning electrical installations and wiring should be strictly complied with, and the same applies to the regulations that may deal with the activities of employees.

It is forbidden to use the neutral conductor, even if earthed, as a protective conductor.

Adequate measures should be taken when more units are installed in parallel.

The control board for A006/l or A006 is designed for an easy operation as the starting and stopping of the motor is automatically controlled by means of the micro-switch located near the operator.

The electrical wiring should therefore include also a low tension line from the surgery to the control board(s) of the Unit(S).

All other Aspirators are to be equipped with a control board that is not normally supplied with them. (Art. Q202)

SECTION B

INSTALLATION, USE AND MAINTENANCE

5 GENERAL

1. All appliance bearing the mark **CE** comply with Directive 93/42/CE (see attached Declaration of Conformity) provided that they are used for the purpose for which they have been manufactured.

2. In order to avoid any danger it is absolutely necessary that all indications in this manual and all other instructions that may be supplied including those contained in the Safety Instructions (SECTION D) be adhered to.

3. Custom's tailored models or variants may differ for technical details. Should any problem arise in this connection, please contact the producer; the model and the Serial Number of the appliance should be given in the enquiry (see them on the plate of the machine) .

6 REQUIREMENTS

6.1 Use

Air ring vacuum pumps are solely to be employed to convey air.

These pumps are designed for continuous use: be careful not to exceed the maximum temperature allowed for the coils when many start and stop operations are carried out in succession or if the ambient temperature is high. Further details in this respect will be supplied on demand. Possible solids or dirt must be removed before they may enter the pump; in this respect please refer to Par 4.2

Whenever possible avoid to install the appliances in the open: please contact 4TEK S.r.l. for any information on this subject.

6.2 Aspiration pressure

If any reduction of diameter in the inlet or outlet piping become necessary, a pressure limiting valve should be inserted.

The total differential pressure is allowed only when the cooling system is in good working conditions. Maximum pressure allowed within the pump is 2 bar absolute. However, if this value is reached, pump performance may be affected.

6.3 Ambient temperature

Max. ambient temperature for these pumps is + 40°C.

7 INSTALLATION

7.1 Transportation

Ensure that the lifting device used to move the appliance is of adequate capacity, referring to the following data:

UNITS WITHOUT CABINET		UNITS WITH CABINET	
A001	10Kg	A004	20Kg
A001/L	8Kg	A004/L	18Kg
A002	14Kg	A005	24Kg
A003	21Kg	A006-OSM	28Kg
A007	30Kg	A006/L-OSM	25Kg

7.2 Assembly

To achieve a smooth and steady running it is necessary that a stable operating conditions be provided. Ensure that cooling of the appliance is not hampered and check that all ventilation slots are unobstructed. An arrow shows the direction of flow. Both Aspirators and Units should never come into contact with inflammable materials (wood, cardboard,...) In Aspirators, when the air can flow freely in and out, the noise level can be reduced by installing a silencer (supplied with special order, Art. S203). A further reduction of the noise level can be obtained by avoiding that the appliance is set on good sound conducting structures (e.g. thin walls, metal sheet panels). Please also ensure that all piping is so laid as to avoid any stress to the appliance.

7.3 Electrical connection

Before performing any work on the appliance, electric power must be cut. Tension and frequency of the mains should match those as indicated on the machine plate. Deviations from these values of the order of +/- 5% of tension and/or +/- 2% of frequency are admitted as they may not cause reduction in power. Earth conductors must be connected to the proper terminal. The motor thermal protection should be adjusted to the current value shown on the plate. Converters may not be used to supply power to the appliances as high frequency tension and current harmonics can cause electromagnetic disturbances and may damage the electronic circuitry. This equipment has been tested and found to comply with the limits for medical devices in IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, the user is encourage to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult that manufacturer or service technician for help

7.4 Starting

Preventive measures should be adopted because temperature of the exposed surface may reach 80°C, this

may cause scalds or burns .
It is furthermore necessary that no temperature sensitive object, like cables or electronic components, should be too close to the appliance.

Nominal rating current in the motor is valid for an ambient and air temperature, at the inlet, of +40°C. If this temperature is +25°C, a 6% current rating is admissible.

8 MAINTENANCE

The Aspirator require a very limited maintenance. Neither the motor, nor the blower require any lubrication. It's however recommended that following operations be carried out.

8.1 Maintenance by the user

External cleaning: the outside of the appliance should be kept free from dust and dirt deposits.

8.2 Maintenance by qualified personnel

Internal cleaning: whenever this should become necessary, please contact 4TEK or send the appliance back to factory for a complete overhaul.

Lubrication: under normal conditions, such as max. ambient and air temperature at the inlet of +40°C and total admitted differential pressure maintained, the following should apply:

1. After approximately 20.000 hours of running or after two and a half years (whichever is the earliest), old grease should be removed from ball bearings and their seats, and new grease should be applied by filling with it about 50% of the space within the bearing and 65% of the seat.
2. Sealed ball bearings must be replaced and their seats must not be greased.

Do not mix different types of grease.

Further information may be obtained contacting 4TEK.

9 TROUBLESHOOTING

Refer to table 2.

10 GUARANTEE

All appliances and parts manufactured by 4TEK are guaranteed free from defects in materials and workmanship for a period of two years.

The conditions of this guarantee are listed in the card that comes with the appliance.

If any assistance is required please contact the local dealer, the representative or directly 4TEK.

WARNING!

- damage caused by faulty electrical or hydraulic connections is excluded from the guarantee;
- the installation of the appliance must be carried out and tested by qualified personnel;
- electrical diagrams supplied with the appliance must be adhered to when connecting the same;
- the characteristics of the power supply should be ascertained beforehand;

- it is mandatory to observe all rules and regulations concerning electrical apparatus and electrical connection with particular reference to the content of par 4.4 of this Manual.

SECTION C SPECIAL INSTRUCTIONS FOR A006L - A006 /OSML - OSM

11 GENERAL FEATURES

A006/L - A006 – OSM/L - OSM is a small centralised unit designed for one surgery.

It consists of an air ring aspirator, a separator, 3 small probes who control the drainage pump and an electronic control board.

12 INSTALLATION

The Unit is best located very close to dental unit, possibly to the left of the chair.

The Unit is manufactured according to the CE standards – to which also the installation shall conform – and can be directly plugged in to the mains as it is equipped with a double fuse switch.

However, before plugging it in, the electrical diagram of the control board as well as the plate data should be examined carefully for their consistency with the mains'; the power supply circuit should include a high sensitivity differential switch.

The tip support tray should be higher than the aspiration inlet of the Unit in order to avoid that secretions may be trapped in the pipe.

The waste outlet, located in the lower part of the Unit, must be connected to the waste through a pipe of the same diameter and should be provided with a siphon. It is very important to check that the drainage line is not affected by other appliances in the dental unit and has an adequate slope.

The air sucked by the Unit is expelled via the holes of the silencer, situated in his lower area.

13 OPERATION

Once the installation is completed, the main switch on the power supply and the Unit switch on the cabinet should be set to ON position. The Unit will start as soon either tip is disengaged from its holder.


Secretions sucked in by the tip **C** (see Pic.2) are separated in the separator **VS** while the air sucked by the aspirator **A** is expelled through the silencer.

Secretions that have accumulated in the separator, are drained through **L** by the drainage pump **P** controlled by the probes, when the canister fills up to 2/3 of its capacity.

14 MAINTENANCE

For all maintenance and cleaning operation, use expendable gloves to protect your hands.

Periodically clean the whole Unit, the tip holder, the aspiration and waste piping.

Pic. 2

Because of the air turbulence induced by the aspiration, blood sucked during surgery generates foams that might interfere with the proper operation of probes, and therefore stop the Unit: adequate anti-foam agents should be used. Always clean the Unit by sucking hot water after each surgical operation.

At the end of each day rinse all parts of the aspiration circuit by sucking liquid agents that have detergent, disinfecting and fungicide action.

It is recommended to wash and disinfect the equipment by immersing firstly the larger tip and then the smaller one in the detergent solution – prepared by dissolving the agent in hot water – and sucking the liquid until a reduction in the aspiration rate is noticed. During this operation, when a certain amount of liquid has been sucked, it is advisable to rise each tip upwards, to improve both aspiration and discharge of the liquid. Repeat this operation two or three times.

At least once every ten days it is necessary to accurately clean the secretion vessel (canister), as well as the valve located in its bottom (viton), and the three probe's terminals. The complex contains the separator vessel can be extracted by pulling out the sleigh by means of the hole in the plate vase-holder.

When it is scheduled that the equipment will be idle for several days, it is necessary to operate it with open tips, in dry conditions – without sucking any liquid at all – for at least 15 minutes before the idle period begins.

This operation serves the purpose to drying up completely the inner parts of the aspirator in order to avoid those typical encrustations, promoted by moisture, that might cause the blocking of the impeller and the stopping of the motor, should their thickness become substantial.

TECHNICAL DATA	A006 -OSM	A006/L- OSM/L
Power	0,37kW	0,415kW
Tension	230V	230V
Frequency available	50Hz 60Hz	50Hz 60Hz
Speed	2800 rpm	2800 rpm
Fuse	CT 5A	CT 5A
Max Vacuum rate	127mbar	110mbar
Suction Capacity	1300lt/min	1250lt/min
Dimensions	40 x 35 x 45h cm	40 x 35 x 45h cm
Weight	28Kgs	25Kgs

15 CONTROL BOARD

A line cord 3.5 m long complete with power plug is supplied with the Unit. The cable is 3x1, Class HO5VV-F, self-extinguishing and includes the protection conduction for earth.

OSM/L – OSM is exactly the same but is not connected with waste pipe (manual drain) and it has accessories: we kindly recommend you to drain the canister 1 every week.

SECTION D SAFETY INSTRUCTIONS

16 INTRODUCTION

16.1 Foreword

This document covers the installation, operation and maintenance of dental aspiration units. Such units are of compact construction and all their potentially dangerous components – such as moving parts or areas exposed to heating – are thoroughly protected.

Even if for easy consultation, following Safety Instructions have not been written to cater for all possible occurrences, they are nevertheless to be adhered to in order to comply with current regulations.

If you wish to receive further information or should solve particular problems not illustrated in this document, please contact 4TEK.

16.2 Definitions

Installation: includes designing and executing the installation of the Unit and its accessories in the premises where it is to be employed as well as its starting.

Operation: includes turning up and using the Unit.

Maintenance: includes checking and preventive steps; maintenance; searching for causes of breakdowns, and fixing of same.

Repairs: includes actions that involve disassembly; replacement of spare parts; reassembling the Unit.


TABLE 2 – TROUBLESHOOTING

TROUBLE	REASON	SOLUTION
1A – The engine does not start, no working noise	1B – In at least two conductors no current flows	1C – Eliminate the problem caused by fuses, terminals or the power supply line
2A – The engine does not start, hum.	2B – In one conductor no current flows. <i>The impeller is blocked.</i> The impeller is faulty. <i>The engine bearing is faulty.</i>	2C – See 1A <i>Open the cover, remove foreign matters, clean.</i> Replace impeller. <i>Replace engine bearing.</i>
3A – The safety switch of the motor disengages after insertion, the current consumption is too high.	3B – Short-circuit in the winding. <i>The engine is overload.</i>	3C – Make a check of the winding. <i>Reduce the operating pressure, if necessary clean the filter and the connecting pipes.</i> See 2A
4A – No vacuum, or insufficient vacuum is generated.	4B – System not waterproof. <i>Wrong direction of rotation.</i> Wrong frequency. <i>The involute profile of impeller's blade is modified because of dirt.</i>	4C – Seal the system. <i>Change the direction of rotation, following the pattern contained in the terminal cover.</i> Rectify the frequency. <i>Clean the impeller or replace it if too worn.</i>
5A – Abnormal grinding noises.	5B – The flow rate is too high. <i>The silencer is dirty.</i> The bearing is without grease or faulty.	5C – Increase the diameter of the hose and clean. <i>Clean or replace silencer's cartridge.</i> Lubricate or replace bearing.

16.3 Symbols

The symbols on the side are used in these Instructions with the meanings as given.

Means that particular attention must be paid to the technical aspects since they may not be immediately obvious even to qualified personnel.

In order to avoid problems that may cause injures persons or damage to property, it is also important that all transportation's, assembly's, operation's and maintenance's instructions be complied with, and that all technical data in manuals or on the same Unit be properly taken into account.

16.4 Qualified personnel

Is the personnel that thanks to its training, experience, knowledge of pertinent specifications, accident prevention regulations and local conditions is authorized to carry out activities such as installation, starting and maintenance of the Units and capable of recognizing and avoiding possible dangers. – A knowledge of First-aid basics is advisable.

17 GENERAL

17.1 Safety notes

Because of their electrical and mechanical properties, if the Units are not used and maintained in accordance with the requirements, or if subjected to inappropriate interventions, they can cause serious damage to people or things.

For this reason it is an essential prerequisite that both design and execution of the installation of electrical and mechanical elements, as well as maintenance, are (with the exceptions described in par 20) carried out and monitored by qualified personnel, who should be fully familiar with these Instructions.

During operation, some parts of the Units are inevitably subjected to dangerous electric current or mechanical stresses. The safe operation of each unit requires of necessity that its transportation and handling operations are adequately carried out and that the Unit is used and

maintained properly. Attention should be paid also to notices and indications on the Unit itself and its components.

17.2 Use

The Aspiration Units under these Instructions are part of groups of appliances employed mainly in the medical sector and their manufacture complies with data shown on their plates and with standards indicated in certificates. This is the reason why these Instructions illustrate only steps to be taken in a normal environment. However, any local regulation relevant to the specific sector should also be adhered to by installing, for instance, special components made to order, or particular additional protections, etc.

18 INSTALLATION

18.1 General safety rules

Completely safe installation and maintenance of the Units are only possible if, with the exceptions under par 20, carried out by qualified personnel who follow the instructions in this document supplied with the Unit and the data indicated in the Unit's plate or its components.

In particular the general standards for installation and safety are to be followed for work on high power appliances as are standards covering the correct use of lifting devices, of tools and protective equipment for the personnel (safety goggles and the like)

Do not reach into the machine through air intakes or discharge ports as the rotating impeller may cause injuries. It is also to be remembered that because of its rotating mass, the machine may continue to turn for several minutes after being shut off.

The impeller may cause injuries even when the machine is switched off, if it is set in motion manually.

18.2 Installation

Under unfavourable operating conditions, parts of the housing and piping may become very hot; depending on the type of installation it may be necessary provide protections against accidental contacts.

It is also to be noted that the fluid being handled may be ejected, at such temperatures, through the discharge ports or the pressure control valves. These openings should be so sited that persons or flammable or explosive materials are not invested by high temperature jets. Temperature sensitive parts such as cables or electronic components are not to be placed next to parts that may become very hot.

The vacuum pumps contained in the Units are suitable only for the transport of air free of dust, or liquids (usually water), non-flammable, non-aggressive and non-explosive.

Solids and impurities are to be intercepted before they enter the Unit by means of filters.

Air ring aspiration Units are air cooled, therefore intake or discharge air must freely circulate. Hot air must not be re-circulated.

No water should penetrate the motor in any circumstance.

Care has to be taken that non accidental contact with moving parts should be ever occur.

18.3 Electrical connection

All wiring must be adequate to the tension and current ratings shown. Connect the Unit according to the diagram in the control box or to that supplied separately for units requiring no such box. In order to avoid danger, connection to the terminals in the box must be carried out professionally.

This implies in particular that:

- the inside of the box is clean and free from remains of old cables;
- the protective or the earth conductor is properly connected;
- all terminals are tight;
- the minimum clearance for discharge in the air is adhered to (pay attention to protruding cable ends);
- unused cables inlet are closed and cover plugs properly screwed in;
- all seals of the box are in a proper condition to maintain the type of protection specified;
- metal surface are clean and lightly lubricated.

CAUTION: Please not change the electrical connector. If you need to change it, we kindly inform you that this has to be done by producer or authorized technician. Changing the power connector with a non-authorized cable could cause harmful interference to other devices in the vicinity and/or malfunctions

18.4 Water and waste lines

Regarding liquid ring Units, before being started, should be connected to the pipes provided so that electrical components cannot be reached by any kind of liquid.

All piping, tanks, valves, taps shall be adequate to the anticipated pressure and temperature and shall be adequate to the type of fluid to be conveyed.

For air exhaust, arrangements contained in Para 4 of this Manual, are to be adhered to whenever the processed air should contain anaesthetics.

18.5 Starting

If the Unit is started up with no connection to the piping or no fastening (for example for test purposes), it may move abruptly or topple over due to the torque reaction of the motor.

For the safe operation of the Unit it is necessary that following conditions are met:

- installation and operating conditions correspond to data on the plate and in the documentation for tension, current, connections, model, type of protection
- the Unit is properly installed and connected to rigid or flexible pipes as the case may be;
- the highness of the location where the Unit is to be installed has been taken into consideration when adjusting the pressure control valves;
- the impeller can rotate without touching anything;
- the direction of rotation is as specified;
- all fastening screws, bolts, fasteners and electrical connections are properly set and made;
- all earth connections and, if any, equipotential connections have been made properly;
- all supplementary devices, if any (temperature probe in coil, anticondensation heater, etc.) are properly connected and operating;
- all measures have been taken to protect against contact with wiring and/or energized parts;
- supplementary fans, if any, are ready for operation and connected for the specified direction of rotation and do not impair the operation of the Unit.

The intake and outlet ports are so located that no foreign material may be sucked in nor be expelled through the discharge port (hazard for eyes and skin, danger of poisoning).

When air is sucked from the atmosphere, the intake ports must be covered by protective devices (gratings and the like) in order to prevent foreign material, including parts of the body or clothing from being sucked in.

To ensure that the safe conditions of the Unit are permanent, the following precautions are recommended at the time of starting and after 500 hours of running :

- check that all screw connections are tight;
- make sure that cables and insulation parts – where accessible – are in good conditions and not discoloured;
- check for mechanical gear noise and vibration on the supports, covers and the elements of the cabinet;
- switch off the Unit if it is not running smoothly or is making strange noises, and start repairs immediately;
- if the Unit is running satisfactorily, check the values for voltage, current and power;
- monitor, as far as possible, the temperature of supports, bearings, etc., until reaching a constant value.

19 OPERATION

19.1 General safety rules

Covers that prevents contact with active or rotating parts or that are required to direct the flow of air for an effective cooling, are not to be opened while the Unit is running.

The noise level of these air ring or liquid ring Units is much below 85 dB (A). Should this level be exceeded beyond certain limits, the Unit must be stopped immediately, a proper inspection has to be carried out followed by suitable corrective measures.

The measures illustrated in Par 18.5 should likewise be taken after protracted idle periods.

20 MAINTENANCE AND REPAIRS

Careful and regular inspections are required to detect any malfunctions at an early stage, before major damage may ensue.

20.1 User's maintenance

Some maintenance may be carried out (see Par 8). While this is done, the user must adhere to these Safety Instructions.

20.2 General safety rules

Before any work is performed on the Unit or its accessories, especially before covers of active parts are removed, the Unit and its accessories are to be properly disconnected from the power source.

Apart from the main electrical circuits, particular attention must be paid to any other auxiliary circuit (see Par 18.2).

Standard safety rules are the following:

- cut the power from source;
- block its reinsertion;
- confirm that no power is available;

protect adjacent parts still energized or provide a barrier.

The above measures are not to be reversed until the Unit is completely reassembled and its maintenance finished

20.3 Dismantling

After screws and bolts are removed, some parts are just held firm by centring fits.

Even if proper dismantling is carried out it is still possible that some heavy parts may therefore suddenly become loose and drop off. This may cause injuries and damage. Take suitable measures to secure all parts being worked on.

If any uncertainty exists we recommend that you contact the producer stating the Unit type and serial number shown on the plate.

20.4 Reassembling

Should prescriptions imposed by the degree of IP protection require that sealants are to be used, joints must be resealed with a suitable non hardening sealant, free from silicone. If gaskets and sealing elements are installed to ensure the degree of protection they are to be inspected and replaced if they are no longer effective.

Tightening torques may be specified in the terminal box for screwed connection of electrical terminals. If these are not complied with, danger may ensue if some cables become loose.

Table 4 - SYMBOLS USED ON THE EQUIPMENT


	CE marked		Production date
	ON		OFF
	Serial number		Read manual
 ATTENTION TURN-OFF THE POWER BEFORE OPENING THE DOOR	Attention: electric shock		

Table 3 - Technical Data

SEE LAST PAGE

4TEK SRL
Viale Dell'industria 66
20037 Paderno Dugnano (MI) ITALY
TEL +39 02 91082211
Fax +39 02 91084668
stefano@4-tek.it gabriella@4-tek.it riccardo@4-tek.it
www.4-tek.it

AIR RING ASPIRATOR UNITS

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ASPIRATORS (without cabinet)
ASPIRATORI (scarenati)
A001/L - A001
A002
A003
A007

UNITS (with cabinet)
UNITA' (caremate)
A004/L - A004
A005
A006/L - A006
OSM/L - OSM

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