CE

ASPIRATORE AD ANELLO UMIDO
WET SUCTION
ASPIRATEUR ANNEAU HUMIDE

RAIN1 – RAIN1C – RAIN1CT - RAIN2

Manuale di installazione ed uso
Installation & User Manual
Notice d'installation et utilisation

4tek S.r.l
Viale Dell'Industria 66
20037 Paderno Dugnano (MI) ITALY
Tel +39 02 91082211 Fax +39 02 91084668
P.IVA C.F. 05133430966
www.4-tek.it
contact@4-tek.it

Sommario

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www.4-tek.it
contact@4-tek.it

Index

1	Instruction warning and caution: please read before operation	21
2	Safety instruction.....	22
2.1	Caution: to prevent injury	22
2.2	Caution: to reduce risk of electrical shock	22
2.3	Danger: to reduce risk of explosion or fire	22
2.4	Description and classification	22
3	Installation notes	23
3.1	Environmental conditions.....	23
3.2	Introduction.....	23
3.3	Transport & storage conditions	23
3.4	Unpacking.....	24
3.5	Mounting.....	24
3.6	Installation location	24
3.7	Suction inlet hose connection	25
3.8	Waste outlet hose connection.....	25
3.9	Exhaust air hose connection.....	25
3.10	Electrical connection.....	26
3.11	Electrical shock hazard.....	26
4	Operation	26
5	Product information	27
5.1	Functional description.....	27
5.2	Wiring diagram.....	30
6	Maintenance.....	30
6.1	Cleaning suction inlet filter.....	30
6.2	Cleaning suction unit	31
7	Disposal	31
8	Trouble shooting chart.....	32
9	Warranty.....	33

1 Instruction warning and caution: please read before operation
 While reading your manual, please pay close attention to areas labeled

WARNING E CAUTION

The description of each is found below

WARNING

Warnings are given where failure to observe instruction could result in injury to people

CAUTION

Cautions are found where failure to observe the instruction could result in damage to the equipment, associated equipment and process

The following symbols of warning will be found on the pump:

Caution / Warning – Refer to accompanying documents

Electrical shock – Attached label at the products

Hot surface – Attached label at the products

	Electric shock		Hot surface
	Alternating current		CE mark
	Attention -consulting manual		Protective earth (ground)
	ON: power		OFF: power
	Date of manufacture		Manufacturer
	Serial Number		Refer to instruction manual
	Dispose of Electric and Electronic Equipment correctly		

2 Safety instruction

2.1 *Caution: to prevent injury*

- Never operate this product if it has a damaged cord or plug. If it is not working properly, has been dropped, damaged or has fallen into water, please contact appropriate service center for examination and repair
- Keep the cord away from heated surfaces. All electrical products generate heat. To avoid serious burns, NEVER touch suction motor during or immediately after operation. The temperature of suction motor can be reach 120°C.
- Never block any air exhaust valve or place it on a soft surface where the openings may be blocked. The air exhaust valve is for ventilation of the motor inside the housing. Keep the air exhaust valve free of lints, dirt and other foreign objects

Blocking air exhaust valve can cause the product to overheat and fire

- Never drop or insert fingers or any other object into any openings
- Use only in well ventilated areas. The motors on all products are totally enclosed fan cooled
- The amalgam separator is guaranteed for 18 months from installation date: after this date it must be replaced.

WARNING

Do not operate the pumps in an atmosphere containing flammable or explosive gases/vapors

- Be sure to properly identify intake and discharge before using
- Protect unit from contaminants and moisture

WARNING

Remove plug from air exhaust port before using

2.2 *Caution: to reduce risk of electrical shock*

- Do not disassemble. Disassembly or attempted repairs if accomplished incorrectly can create electrical shock hazard. Refer servicing to qualified service agencies only
- Do not use this product in or near area where it can fall or be pulled into water or other liquids
- Do not reach for this product if it has fallen into liquid. Unplug immediately
- Never operate this product outdoors in the rain or in a wet area

2.3 *Danger: to reduce risk of explosion or fire*

- Do not use this unit in or near explosive atmospheres or where aerosol (spray) products are being used
- Do not use this product near flames

2.4 *Description and classification*

The equipment is a wet suction used during dental treatment. It consists of an air ring aspirator, a separator and a control board; for model RAIN1CT the amalgam separator is

included.

When the dentist put the cannulas up, the air ring suction starts to run so it sucks secretions, water and blood from the patient's mouth. They will go inside of the separator that, being assembled on the same tree of the electrical motor, will use the centrifugal force to separate water and liquids from air: so water will be discharged in waste line and air will be expelled outside of the building through the exhausting pipe.

In the case of the RAIN1CT, an amalgam separator compliant with ISO 11143:2018 is interposed between the discharge line and the device: this separator cannot be opened but it must be removed after 18 months from installation and replaced with a new one. The removed one must be handed over to a company authorized to collect and discharge amalgam and special waste.

- Class I equipment according to EU REG 745/2017
- Ordinary equipment (IPX0)
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide
- Continuous operation

This equipment has been tested and found to comply with the limits for medical devices in IEC 60601-1-2. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encourage to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device
- Increase the separation between the equipment
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected
- Consult that manufacturer or field service technician for help

3 Installation notes

3.1 *Environmental conditions*

The pump is rated for indoor use only.

Maximum altitude 2000 meters.

Operating temperature range 10°C to 40°C.

Maximum relative humidity of 80% for temperature up to 31°C decreasing to 50% at 40°C.

Rated for +/- 10% of supply voltage.

Pollution degree 2.

Installation category II

3.2 *Introduction*

This manual has been compiled not only for the care and maintenance of the suction system now in your possession, but as a helpful reference and guide to prevent many problems which can occur if used improperly.

3.3 *Transport & storage conditions*

Temperature: -20°C to 40°C

Humidity: 0% to 90%

Pressure: 700hPa to 1060hPa

3.4 Unpacking

Carefully remove the machine from the shipping case. Preserve all paperwork for future reference. If damage has occurred from shipment a claim must be filed with the carrier immediately; preserve the shipping carton for inspection by the carrier. If you are required to communicate with your dealer, be sure to include serial numbers for quick identification.

3.5 Mounting

Rubber feet are attached to the product. Rubber feet are excellent for applications involving a semi-flexible surface such as a bench top; they help to isolate noise and eliminate creeping. However, the horizontal position of the motor shaft should be preferred.

3.6 Installation location

The suction system should be located preferably in a clean, dry and well ventilated area. Please be sure not to block the ventilation ports located on the motor housing. The unit should be placed where the surrounding temperature remains between 10°C and 40°C (50°F and 104°F). Always check to insure the location chosen is protected from direct or indirect moisture contact.

Place suction outlet valve higher position than drainage to prevent backflow problem.

ATTENTION: The maximum allowed distance between the suction and the dental chair is 6m for RAIN1, RAIN1C, RAIN1CT and 8m for RAIN2.

Hose connection


Diagram of Suction tubing of RAIN1/RAIN1C/RAIN1CT


Diagram of Suction Tubing of RAIN2

Please be sure not to be bended any of suction inlet hose and outlet hose. This can causes loosing vacuum power or damage of suction motor. Also use exhaust air hose with 40mm diameter hose.

Exhaust air temperature increases significantly and discharge drastically too hot for most plastic piping. Therefore, metal piping or heat-resistant pipe is recommended for at least the first five (5) to eight(8) feet (1.5 ~ 2.5 meters) from the air exhaust valve on the discharge side. In addition, this piping **MUST** be guarded and marked “**DANGER–HOT–DO NOT TOUCH.**”


CAUTION

Do not use hoses such as rubber or PVC hose which are not resistant to chemicals or disinfectants or insufficient flexibility. Only use PVC flexible spiral hoses with integrated spiral or hoses of equivalent material

3.7 Suction inlet hose connection

- Remove the cover of suction inlet and push-in tubing connector in the inlet hole
- Connect suction inlet hose (Ø40mm spiral hose) at tubing connector and tighten with teflon insulation tape

3.8 Waste outlet hose connection


- Always use spiral hose with Ø20mm interior
- Connect hose to the outlet hole with clockwise
- Please be tighten with teflon
- Always be higher position than drainage
- For RAIN1CT connect the outlet hose of the device to the entrance of the amalgam separator and the outlet hose of the amalgam separator to the discharged pipe

3.9 Exhaust air hose connection

- Turn left exhaust air outlet hole clamp
- Connect exhaust air outlet hose to the hole (strongly push-in)
- Turn right exhaust air outlet hole clamp for fasten
- Always place end of the outlet hose out of building or window to prevent noisy

problem, hot temperature and bad odor

3.10 Electrical connection

The mains connection must only be made by well-trained electricians.

To insure safe installation, check local current and voltage. Serious damage may occur to the motor if it is connected improper voltage.

CAUTION: Please not change the electrical connector. If you need to change it, we kindly inform you that this has to be done by producer or authorized technician. Changing the power connector with a non-authorized cable could cause harmful interference to other devices in the vicinity and/or malfunctions.

3.11 Electrical shock hazard

Disconnect electrical power at the circuit breaker or fuse box before installing this product. Install the system in a location where it will not come in to contact with water or other liquids. Install the system in a location protected from the weather. Electrically ground the system. Failure to follow these instructions can resulting death, fire or electrical shock.

4 Operation

Check list before start up the system

- Check if there is any foreign material in the cooling fan part of the motor
- From the motor side of the system, verify if the blower is rotating in the direction of clockwise. On the system powered by a 3-phase electric power, change the connection of any two (2) wires to reverse blower rotation
- Check any abnormal noise or vibration of the suction system. If this happen, refer to the trouble shooting guide on this manual. After removing causes of the noise and vibration, re-start the system.

5 Product information

5.1 Functional description


RAIN1 / RAIN1C/RAIN1CT

1. Vacuum connection tube
2. Air/liquid separator
3. Motor impeller chamber
4. Separation unit case
5. Waste water outlet
6. Auto by-pass drain tubing
7. Rubber foot
8. Suction inlet hose
9. Filter
10. Filter housing
11. Electric control unit
12. Power switch
13. Cooling fan cover
14. Exhaust air outlet

Amalgam separator according to ISO11143:2018


RAIN 2

1. Cooling fan cover
2. Electrical box of motor
3. Motor
4. Impeller chamber
5. Electrical box
6. Exhaust air outlet
7. Waste water outlet hole
8. Auto by-pass draining
9. Suction inlet hose
10. Filter
11. Filter housing
12. Power cord
13. Electrical switch with thermic

TECHNICAL DATA SHEETS		
	RAIN1/RAIN1C/RAIN1CT	RAIN2
Main voltage (VAC)	230 ~ 1 phase	230 ~ 1 phase
Frequency (Hz)	50 / 60	50 / 60
Power (Kw)	0.55 / 0.63	0.7 / 0.83
Vacuum level (mBar)	170 / 200	220 / 290
Suction motor rotation speed	2800 / 3400	2800 / 3400
N. of chairs	1	2
Weight (Kg)	RAIN1 19 RAIN1C 20 RAIN1CT 21	20
Dimensions	RAIN1 28 x 50 x 55h cm RAIN1C 33 x 57 x 60h cm RAIN1CT 28 x 50 x 78h cm	34 x 43 x 54h
Sound level (dB)	RAIN1 53 RAIN1C 50 RAIN1CT 53	55
DIMENSIONS OF HOSE CONNECTIONS: <ul style="list-style-type: none"> - Suction: 40mm (internal) - Draining: 20mm (internal) - Exhaust air: 40mm (internal) - For Amalgam separator: D.in: 25mm D.out: 32mm 		

5.2 Wiring diagram

6 Maintenance

6.1 Cleaning suction inlet filter

- Pull out filter part from machine in counter-clockwise direction of cover.
- Clean filter in running water for removing particle or debris.

(Recommending it cleaned at least twice per month for avoiding reduction of suction power)

- Re-connect filter

For RAIN1CT unit with amalgam separator RCT1, since it cannot be opened, it does not require any maintenance but the separator must necessarily be replaced 18 months after installation.

CAUTION

If the machine installed into the used suction tubing line, old debris and dirty particle may clogged filter. In case of this, please clean filter frequently for some period of time until tubing line is cleaned.

6.2 Cleaning suction unit

After every treatment: For avoiding bad odor and reducing infection possibility, clinic should suck 1 glass of water before the next treatment

After finishing last treatment of the day – For avoiding bad odor, reducing infections possibility and assured a properly working of the unit, even to obtain extension of the warranty, dentist must clean suction holes sucking **NON-FOAMING** and **DISINFECTANT** liquids.

CAUTION

DO NOT use a foaming cleaning
It causes backflow of water to the suction motor and can damage it

7 Disposal

CAUTION

This product may be contaminated after using years. Please make the disposal organization aware of this in order that they can take the proper safety actions

Uncontaminated plastic components of the product may be recycled

The built-in electric circuit board and other electric components should be disposed as electric waste. Other metal components may be disposed as metallic waste. If the product

is returned to the local dealer or 4TEK S.r.l., all connections should be closed so that they are water-tight.

For RAIN1CT devices equipped with amalgam separator, it must be replaced with a new one 18 months after installation: once the old separator has been disconnected from the system, after closing the inlet and outlet unions with the appropriate caps supplied with the separator, it must be sent to a company authorized to collect and discharge the amalgam.

8 Trouble shooting chart

This trouble-shooting is carried out by qualified technicians only.

Problem	Possible reasons	Remedy
Unit does not work	Mains supply is not connected Incorrect electrical connection or power source	Check mains power, fuses and reset circuit breaker if necessary
Unusual noise increased	Solid particle or debris is coming into the suction motor	Dismantle unit and clean the impeller of the motor
Excessive vibration	Damaged impeller Motor and/or impeller are dirty	Dismantle unit and clean the impeller of the motor
Unit is very hot	Wrong wiring Low voltage Suction inlet is clogged Operating pressure or vacuum is too high	Check wiring Supply proper voltage Clean primary filter Install a relief valve and pressure or vacuum gauge
Suction is too low	Suction inlet is clogged Leaking problem in the suction plumbing	Clean primary filter Check all suction tubes, hoses and connections. Replace if necessary
Motor overload	Low voltage	Check power source Check wire size and wire connections

9 Warranty

4TEK S.r.l Suction System, when properly installed and operated under normal conditions of use, are warranted by 4TEK S.r.l. to be free from defects in material and workmanship for a period of 24 months from the date of receiving the 4TEK's guarantee card or however no more than 36 months after the date of 4TEK invoice.

In order to obtain performance under this warranty, the buyer must promptly (no later than 30 days after discovery of the defect) give written notice of the defect to 4TEK S.r.l. (Fax: +39 02 91084668 or E-mail: contact@4-tek.it) or to a representative in his country. Customer is responsible for freight charges to 4TEK S.r.l. in all cases.

This warranty does not extend to any goods or parts which have been subjected to misuse, lack of maintenance, neglect, damage by accident or transit damage.

4TEK S.r.l. is not responsible or liable for indirect or consequential damages of any kind however, including but not limited to those for use of any products, loss of time, inconvenience, lost profit, labor charges, or other incidental or consequential damages with respect to persons, business, or property, whether as a result of breach of warranty, negligence or otherwise.

This warranty can be modified only by authorized 4TEK personnel, by signing a specific written description of any modifications.

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Seule la société 4tek S.r.l. peut modifier cette garantie.