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
English


Table of Contents

1	General.....	4
1.1	Contact information	5
1.2	Structure of the document.....	5
1.2.1	Labeling of information.....	5
1.2.2	Formats and symbols.....	5
1.2.3	Abbreviations/codes used	6
1.3	Trade names / registered trade marks of the manufacturers	6
2	Safety notices.....	7
3	Technical description.....	9
3.1	Task.....	9
3.2	Functionality	9
3.3	Structure of Midwest T1 / T2 Energo.....	9
3.4	Technical data	10
3.5	Product labeling.....	11
3.6	Standards and approvals	11
4	Preparation.....	12
4.1	Initial start-up and longer breaks in use	12
4.2	Prior to starting the work day.....	12
4.3	Before each patient	12
5	Operation.....	13
5.1	Replacing the instrument.....	13
5.2	Inserting and removing burs.....	14
6	Postprocessing.....	15
6.1	After each treatment session.....	15
6.2	At the end of the work day.....	15
7	Reprocessing	16
7.1	Manual cleaning with a brush	16
7.2	Care.....	17
7.2.1	Lubricate mechanical parts automatically	17
7.2.2	Lubricate mechanical parts manually	17
7.2.3	Lubricate the push button chuck	18
7.2.4	Clean cooling spray nozzle opening	18
7.2.5	Cleaning the fiber optic surfaces.....	18
7.3	Sterilization.....	19

8 **Maintenance**..... 20

8.1 Check the FG clamping system with a chuck tester..... 20


8.2 Check the FG clamping system without a chuck tester..... 20

9 **Spare parts and consumables**..... 21

10 **Disposal**..... 22

1 General

Scope of the instructions for use

These instructions for use are valid for the following instruments:

- Midwest® T1 Energo S; S = short, short instrument (A)
- Midwest® T1 Energo (B)
- Midwest® T2 Energo (B)
- Midwest® T2 Energo M; M = Mini, mini head (C)

Any mention of Midwest® T1 / T2 Energo in the following text refers to all above-mentioned products.

Before you begin....

⚠ CAUTION! Risk of injury!

1. Read the instructions for use prior to using Midwest® T1 / T2 Energo.
2. Use Midwest® T1 / T2 Energo only for the applications described in the instructions for use.
3. Observe the applicable hygiene standards, occupational safety regulations, and accident prevention measures for Midwest® T1 / T2 Energo in your country.

Store documents

Keep the instructions for use so that they are always at hand in case you or another user requires information at a later point in time. Save the instructions for use on the PC or print them out.

If you sell the unit, make sure that the instructions for use are included with it either as a hard copy or on an electronic storage device so that the new owner can familiarize himself with its functions and the warning and safety notices.

"Download Center" for technical documents

We have set up a "Download Center" for the technical documents at dentsplysirona.com/ifu. From here, you can download these instructions for use along with other documents. Please complete the online form if you would like a hard copy of the instructions for use or operator's manual. We would be happy to send you a printed copy, free of charge.

"Customer Support Portal"

You can obtain additional product information, tutorial videos and other self-help tips at the "Customer Support Portal" at www.dentsplysirona.com/csp.

Help

If you still need assistance, despite carefully reading this technical document and the information in the "Customer Support Portal", contact your dental depot.

Target group

Midwest® T1 / T2 Energo is intended only for use by trained dental personnel in dental practices or laboratories. Midwest® T1 / T2 Energo meet the ISO 14457 standard.

⚠ **CAUTION! Risk of injury!** The instruments may be operated only by trained personnel.

Intended use

Intended for use in dental practices to prepare dental cavities for restorations such as fillings, and for cleaning teeth.

Contraindications

none

Mandatory reporting

Notify the manufacturer and the responsible national authority immediately of any serious incident that has occurred with the product.

1.1 Contact information

Dentsply Sirona Product service

Log in to register your units and make service requests:
<https://dentsplysirona.service-pacemaker.com/>

Manufacturer's address

Sirona Dental Systems GmbH
Fabrikstrasse 31
64625 Bensheim
Germany

Tel.: +49 (0) 6251/16-0
Fax: +49 (0) 6251/16-2591
E-Mail: contact@dentsplysirona.com
www.dentsplysirona.com

1.2 Structure of the document

1.2.1 Labeling of information

Warnings

> To prevent injuries, please observe the warnings.

Warnings are labeled as follows:

⚠ **DANGER!** indicates a danger **leading** to death or serious injury if not avoided.

⚠ **WARNING!** indicates a danger that **may lead** to death or serious injury if not avoided.

⚠ **CAUTION!** indicates a danger that **may lead** to injury if not avoided.

Instructions for use

> To prevent material damage and additional expenses, please observe all instructions for use.

Instructions for use are labeled as follows:

NOTICE! indicates measures for the prevention of material damage.

IMPORTANT: indicates information on the avoidance of additional expenses and other important information.

Tip: indicates information for facilitating work.

1.2.2 Formats and symbols

The formats and symbols used in this document have the following meaning:

✓ Requirement 1. First action step 2. Second action step or > Alternative action ↪ Result > Individual action step	Identifies an action sequence with requirement and result.
Use of formats and symbols [→ 5].	Identifies a reference to another text passage and specifies its page number.
• List	Identifies a list.

1.2.3 Abbreviations/codes used

FG	Friction grip bur instrument
CA	Contra-angle handpiece bur instrument
HP	Straight handpiece bur instrument

1.3 Trade names / registered trade marks of the manufacturers

The designations used in this operating instruction with reference to the instruments are trade names or registered trade names of the manufacturers:

KaVo, Biberach, Germany:

- INTRAmatic
- INTRAmatic LUX

2 Safety notices

Obligations of the user	Protect yourself, the patient, and any third parties from danger. Observe the safety information, the intended use, the applicable hygiene standards, occupational safety regulations, and accident prevention measures.
Preventing the spread of infections and cross contamination	▲ CAUTION! Risk of contamination! Prevent the spread of infection and cross contamination between patients, users and third parties. Before initial start-up and after every treatment, the product and accessories must be reprocessed, see "Reprocessing" section. Take the appropriate hygiene measures, e.g., wear protective gloves.
Filtered water	▲ CAUTION! Risk of burns! Use filtered water only (< 50 µm).
Emitted cooling air	NOTICE! The cooling air emitted by the coupling of the motor must have a flow rate of 1.5 - 10 l/min.
Sufficient cooling	▲ CAUTION! Risk of burns! Insufficient cooling leads to overheating of the preparation site and to damage to the tooth substance. Ensure that the water content is ≥ 50 mL/min.
Instrument head overheating	▲ CAUTION! Risk of burns! If the instrument is defective, the area around the instrument head may heat up. ▲ CAUTION! Risk of burns! A damaged or incorrect bur (bur outside the specification) can cause the area around the instrument head to heat up. Use only burs, as described in section "Technical data [→ 10]", that are undamaged.
Fall damage	NOTICE! Dropping an instrument can lead to fall damage of the instrument. Check the instrument for proper functioning after a fall.
Malfunction or damage	▲ CAUTION! Risk of injury/Risk of burns! In case of malfunction, an unusual or different operating sound, or damage, stop using the instrument immediately and notify your dental depot or the manufacturer.
Repair	▲ CAUTION! Risk of burns! Do not repair the instrument yourself. Safe operation cannot be guaranteed if proper repairs have not been carried out. Contact the dental depot / manufacturer.
Spare parts and accessories	▲ CAUTION! Risk of injury! Use only original Dentsply Sirona parts or parts approved by Dentsply Sirona. Safe operation is not guaranteed for parts that have not been approved by Dentsply Sirona. The warranty does not cover defective products for which original spare parts were not used.
First aid measures in the case of lubricant accidents	General information: Immediately remove any clothing soiled by the product. After inhalation: Supply fresh air; consult doctor in case of complaints. After skin contact: If skin irritation continues, consult a doctor.

After eye contact: Rinse opened eye for several minutes under running water.

After swallowing: If symptoms persist, consult a doctor.

Operating conditions: Temperature: +10 °C – +30 °C.

If you have any questions, please contact your dental depot or the manufacturer.

3 Technical description


3.1 Task

The instrument transmits the driving power and speed of the electric motor or air motor to the preparation instrument.

3.2 Functionality

The transmission ratio of the instrument distorts the operating speed of the preparation instrument.


3.3 Structure of Midwest T1 / T2 Energo

A	Instrument head
B	Push button
C	Cooling spray outlet
D	Opening of the chuck system
E	Light aperture
F	Color coding for gear ratio
G	ISOslider


3.4 Technical data

Contra-angle handpiece

Energy model	Gear ratio			Handpiece length	ISO 3964
	1:5	1:1	5:1		 L = max. 31,8 mm M = max. 24,8 mm S = max. 22,8 mm
Midwest®T1 Energo	X	X	X	long	L, M, S
Midwest®T1 Energo S	X	X	--	short	S
Midwest®T2 Energo	X	X	X	long	L, M, S
Midwest®T2 Energo M	X	--	--	long	L, M, S

Technical data for all Midwest® T1/T2 Energo				
Gear ratio		1:5	1:1	5:1
Color coding for gear ratio		Red	Blue	Green
Maximum motor speed in rpm		40,000	40,000	40,000
Maximum operating speed in rpm		~200,000	~40,000	~8.000
Clamping system		FG	CA	CA
Internal cooling media		X	X	X
Back suction stop		X	X	X
Light function		X	X	X
Instrument coupling		INTRAmatic LUX ISO 3964		

Bur

		FG	CA
	Shank diameter in mm	1.590 - 1.600	2.334 - 2.350
	Maximum total length in mm	25	34
	Maximum working diameter in mm	2.1	-
	Gripping length in mm	≥ 11	≥ 11
	Standard	ISO 1797, Type 3 ISO 2157	ISO 1797, Type 1

⚠ CAUTION! Risk of burns! A damaged or incorrect bur can cause the area around the instrument head to heat up. Use only burs, as described above, that are undamaged.

Treatment center

	Midwest® T1 / T2 Energo
Spray air pressure in bar	2.7
Spray water pressure in bar	2.0
Maximum water temperature in °C	40
Recommended water content in spray in mL/min	> 50


3.5 Product labeling

Sterilizable in a steam sterilizer (autoclave) at the temperature specified

Serial Number

Date of manufacture

Purchase order number

Data matrix code (here: an example)

Content of the data matrix code:

+E2768752050/\$+2002/16D20180829A

(A)

(B)

(C)

(D)

A Manufacturer code (here: E276)

B REF number (here: 875205)

C Serial no. (here: 2002)

D Date of manufacture (YYYYMMDD)

3.6 Standards and approvals

This product bears the CE marking in accordance with the provisions of the Regulation (EU) 2017/745 of April 5, 2017 concerning medical devices.

This product is a medical device.

4 Preparation

4.1 Initial start-up and longer breaks in use

- > Clean and sterilize the instrument and accessories prior to startup [→ 19].
- > Clean and maintain the instrument after longer breaks in use [→ 16].

4.2 Prior to starting the work day

- > **⚠ CAUTION! Risk of contamination!** Purge the water and air channels for 30 seconds.

4.3 Before each patient

1. **⚠ CAUTION!** Purge the water and air channels for 30 seconds.
2. **NOTICE!** Stop the motor.
3. **⚠ CAUTION!** Do not attach the instrument until the motor has stopped. [→ 13]
4. Insert the preparation tool [→ 14].
5. **⚠ CAUTION!** Set a sufficient amount of cooling water (> 50 mL/min).
6. **⚠ CAUTION!** Use filtered water only (< 50 µm).
7. **⚠ CAUTION!** Check the nozzles for blockages and lime deposits, for example, and clean the nozzles if necessary [→ 18].
8. Check the fiber optic surface and clean it if needed [→ 18].

5 Operation

⚠ **CAUTION! Risk of burns!** Insufficient cooling leads to overheating of the preparation site and to damage to the tooth substance. Ensure that the water content is ≥ 50 mL/min.

⚠ **CAUTION! Risk of burns!** Always operate contra-angle handpieces marked with a red dot with cooling spray!

⚠ **CAUTION! Risk of burns!** Never pull the patient's cheek back with the contra-angle handpiece while the motor is running! This would actuate the push button, thus creating a risk of burning the patient's oral mucosa.

⚠ **CAUTION! Risk of burns at push button!** Press the push button only when the instrument has come to a stop.


⚠ **CAUTION! Risk of burns!** For use and operating data of the different preparation instruments, please refer to the manufacturer's specifications.

5.1 Replacing the instrument

NOTICE! The instrument should be attached or removed only when the motor is at standstill.

NOTICE! Do not use instruments with light function on motors without light function.

Attaching the instrument

✓ The motor is at a standstill.

1. Insert the instrument until it snaps into place.
2. When using a motor with a light groove (A): turn the instrument until the ISO slider clicks into place.

Removing the instrument

✓ The motor is at a standstill.

➤ Detach the instrument. Do not pull on the supply hose.

5.2 Inserting and removing burs

⚠ **CAUTION! Risk of burns!** Ensure that the push button is not jammed and thus not pressed in during treatment.


⚠ **CAUTION! Risk of injury!** Use only sharp, undamaged, ISO 1797-compliant preparation tools with a metallic shaft.

⚠ **CAUTION! Risk of injury!** A loose or partially removed preparation instrument can detach itself from the head or break off. There is a risk of injury! Use the instrument only if the preparation instrument is attached at least 10 mm deep or up to the stop and is clamped securely in place.

⚠ **CAUTION! Risk of injury!** Use only perfectly functioning burs that comply with the respective manufacturer's specifications. Also observe the speed recommended by the manufacturer during preparation.

Insert Friction Grip bur (FG)

- ✓ The motor and bur are at a standstill.
- 1. Press the push button and insert the bur as far as it will go.
- 2. ⚠ **CAUTION! Risk of injury!** Pull on the bur to check that it is held firmly in place.


Remove Friction Grip bur (FG)

- ✓ The motor and bur are at a standstill.
- Press the push button and remove the bur with tweezers.
- ⚠ **CAUTION! Risk of injury!** Press the push button chuck only when the instrument is at a standstill.

Insert preparation instrument (CA)

- ✓ The motor is at a standstill.
- 1. Insert the preparation instrument without pressing the push button.
- 2. Turn the preparation instrument back and forth gently until it snaps into place.
- 3. ⚠ **CAUTION! Risk of injury!** Pull and turn the preparation instrument to check that it is firmly seated.

Remove preparation instrument (CA)

- ✓ The motor and preparation instrument are at a standstill.
- Press the push button and remove the preparation instrument.

6 Postprocessing

6.1 After each treatment session

⚠ **CAUTION! Risk of contamination!** Perform reprocessing immediately after treatment, after 1 hour at the latest.

- ✓ The preparation instrument is at a standstill.
- ✓ Wear suitable protective clothing.
- 1. Purge the water and air channels on the treatment center for 30 seconds.
- 2. Remove the preparation instrument with tweezers.
- 3. Detach the instrument from the motor.
- 4. Transport the instrument to the hygiene room in a suitable transport container.
- 5. ⚠ **CAUTION! Risk of contamination!** Conduct hygienic reprocessing [→ 16]. Observe the relevant national/local regulations.

6.2 At the end of the work day

NOTICE! Do not leave any instruments on the motor overnight, in order to prevent oil from leaking into the electric motor. Never lubricate the electric motor.

7 Reprocessing

⚠ **CAUTION! Risk of contamination!** Perform reprocessing immediately after treatment, after 1 hour at the latest.

⚠ **WARNING! Risk of infection!** Reprocess your instrument as described in the following steps.

⚠ **CAUTION!** Do not use agents containing hydrogen peroxide (H₂O₂) or chloride.

7.1 Manual cleaning with a brush

IMPORTANT: When carrying out manual reprocessing, the national/ local regulations must be followed. These must be checked beforehand.

NOTICE! **Never** clean in an ultrasonic bath!

⚠ **CAUTION! Risk of injury! Never** immerse in disinfectants!

⚠ **CAUTION! Risk of contamination!** Use a soft, clean, and disinfected brush for thorough cleaning.

IMPORTANT: The Midwest Handpiece Cleaner is a surfactant.

✓ Wear suitable protective clothing.

✓ The air and water lines are purged for 30 seconds.

1. Thoroughly brush the product under running water (<38 °C, <100 °F, at least drinking water quality) for at least 10 seconds.

2. Use a soft, clean, and disinfected brush for thorough cleaning. Midwest Handpiece Cleaner Plus can be used to support the process.

The instruments must be rinsed to remove disinfectant residue for the subsequent steps.


3. Check at good lighting (min. 500 lux) if the product is clean after reprocessing.

4. If there is any contamination, repeat the process.

5. Lubricate mechanical parts automatically [→ 17] or manually [→ 17].

6. Lubricate the push button chuck [→ 18].

7. Perform sterilization [→ 19] afterwards.

7.2 Care

7.2.1 Lubricate mechanical parts automatically

Interval

- After each patient


Required accessories

- Midwest Automate

Procedure

Lubrication of the instrument can be performed by the Midwest Automate Handpiece Maintenance System.

For further details and instructions, refer to the operating instruction supplied with unit.

If another system is used, the maintenance equipment must be approved by its manufacturer for the lubrication of dental handpieces.

⚠ CAUTION! Risk of burns! Overlubrication of the instrument results in overheating of the instrument head.

7.2.2 Lubricate mechanical parts manually

Interval

- After each patient

Required accessories


- Spray adapter
- O-ring for spray adapter
- Midwest® Plus Spray

⚠ CAUTION! Risk of injury! Use only Midwest® Plus Spray.

⚠ CAUTION! Risk of burns! Overlubrication of the instrument results in overheating of the instrument head.

Procedure

- ✓ The spray adapter is disinfected.
 - ✓ The O-ring on the spray adapter is intact.
1. Fit the spray adapter onto the nozzle of the spray can.
 2. Insert the instrument until it snaps into place and hold it.
 3. Spray into the instrument for 1-2 seconds.
IMPORTANT: Hold the spray can upright.
 4. Wipe away any spray that runs out with a clean, lint-free cloth.
 5. Repeat the process until the spray that runs out is clear.
 6. Purge the instrument with compressed air (3 bar) for 10 seconds.

7.2.3 Lubricate the push button chuck

Use Midwest Plus Spray on the push-button chuck to remove deposits and ensure proper functioning of the clamping system.

⚠ CAUTION! Risk of injury! Make sure you observe the following intervals.



Interval

- At least once a week

Required accessories

- Midwest® Plus Spray

⚠ CAUTION! Risk of injury! Use only Midwest® Plus Spray.

Procedure

- ✓ The instrument is clean.
- 1. Press the contra-angle handpiece head firmly against the spray can nozzle using the chuck.
- 2. While pressing the push button, spray through the chuck for 1-2 seconds.
IMPORTANT: Hold the spray can upright.
- 3. Wipe away any spray that runs out with a clean, lint-free cloth.
- 4. Repeat the process until the spray that runs out is clear.
- 5. Purge the instrument with compressed air (3 bar) for 10 seconds.

7.2.4 Clean cooling spray nozzle opening

⚠ CAUTION! Risk of burns! Regularly clean the openings of the cooling spray nozzles with the cleaning wire.


Interval

- At least once a week

Procedure

- ✓ The instrument is clean and lubricated.
- 1. Carefully clean the nozzle openings by regularly running a cleaning wire through them.
- 2. Let the instrument run briefly with cooling spray.

⚠ CAUTION! Risk of burns! Insufficient cooling leads to overheating of the preparation site and to damage to the tooth substance. Ensure that the water content is ≥ 50 mL/min.


7.2.5 Cleaning the fiber optic surfaces

Interval

- Check and clean if needed

Procedure

- ✓ The instrument is clean and lubricated.
- ✓ A Q-tip/cloth and alcohol are available.

1. In order to avoid scratching the surfaces, blow off any dirt particles with air from a dry syringe.
2. Wipe the surfaces with a Q-tip or a soft cloth and alcohol.

7.3 Sterilization

⚠ CAUTION! Risk of contamination! Make sure you observe the following intervals.

Intervals

- Prior to initial operation
- Prior to each subsequent use

Wrapped sterilization

- ✓ The instrument is clean [→ 16].
- ✓ The instrument is lubricated [→ 17].
- ✓ The instrument is packed in packaging suitable for sterilization and storage (FDA cleared, paper/plastic composite packaging according to ISO 11607).
- Sterilize the instrument in the steam sterilizer with saturated water vapor.

Temperature: 132 °C (270 °F) / 135 °C (275 °F)

Overpressure: 1.87 bar (27.12 psi) / 2.13 bar (30.89 psi)

Packaging	Steam sterilization	Holding time at least	Drying time at least
Wrapped	Gravitation process	15 minutes at 132 °C (270 °F)	15 minutes
Wrapped	Gravitation process	10 minutes at 135 °C (275 °F)	30 minutes
Wrapped	Dynamic air removal	4 minutes at 132 °C (270 °F)	20 minutes
Wrapped	Dynamic air removal	3 minutes at 135 °C (275 °F)	16 minutes

NOTICE! Do not exceed 140 °C (284 °F) even during the drying phase.

After sterilization

1. Remove the instrument from the steam sterilizer.
 - ⚠ CAUTION!** The instrument may still be hot.
 - NOTICE!** Do **not** attempt to accelerate the cooling process by immersing the instrument in cold water. This would damage your instrument.
2. Perform a visual check to see if the packaging is dry after sterilization.
3. Store all instruments so that they are protected from contamination.
4. Sterilize again once the storage period has elapsed.

8 Maintenance

8.1 Check the FG clamping system with a chuck tester

Interval


- at least once a month

Required accessories

- Chuck tester

If the chuck tester is not available, follow the steps described in "Check the FG clamping system without a chuck tester [→ 20]".

Procedure

- ✓ The expiry date of the chuck tester (A) has **not** elapsed.

1. Insert the chuck tester into the FG clamping system.
2. Tighten the chuck tester until the marking ring appears (withdrawal force: 22 N).
3. Check if the chuck tester slips out of the chuck before the marking ring appears.

✗ If no: The chuck is functional.

✗ If yes: **CAUTION! Risk of injury!** The chuck is defective and the secure fit of the preparation tool cannot be guaranteed. Stop using the instrument. Have the instrument's clamping system replaced in a workshop authorized by Midwest.

Tip: Record the time and result of the check for your own information.

8.2 Check the FG clamping system without a chuck tester


Interval

- at least once a month

Required accessories

- Large head diamond bur

Procedure

- ✓ You wear suitable protective gloves.

- ✓ A diamond bur with a large head is ready.

1. Insert a large head diamond bur into the FG clamping system.
2. Carefully pull on the bur (at least 5 lbs) and check to make sure the bur stays in place.
3. If the bur moves or pulls out, do not use and send the instrument in for service.

Tip: Record the time and result of the check for your own information.

9 Spare parts and consumables

	REF		REF
Chuck tester	33 27 793	Midwest Plus Spray	38 00 80M
Cleaning wire for spray nozzles	24 00 232	Midwest Plus Handpiece Cleaner	38 01 40
Spray adapter	89 17 858		
O-ring for spray adapter	70 36 353		

10 Disposal

- According to current information, the product does not contain any substances that are hazardous to the environment.
- ⚠ **CAUTION! Risk of contamination!** Disinfect the device prior to disposal.
- Observe the applicable disposal regulations for your area.

We reserve the right to make any alterations which may be required due to technical improvements.

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© Sirona Dental Systems GmbH
SIRONA Dental Systems GmbH