

True 1-hour
Single Visits with
Lithium Disilicate
Restorations

Amber® Mill Direct


SPEED / DURABILITY / REAL GRADATION

Pre-cyrstallized lithium disilicate block that achieves excellent restorations without the need for an oven.



True Single-Visit Results

We solve the challenges faced with indirect millable restoration materials





YAO-LIN TANG, DDS AestheticLab

"Amber Mill Direct has all the advantages of lithium disilicate ceramics. Its power, however, are the beautiful smooth margins without the need for firing – an invaluable CAD block every dentist should have in their office"



CHRITIAN PETRI, CDT Oral Design Clinic

"No glaze, no stain, just MILL & POLISH, and the final restoration is ready. Anybody can do it, so don't wait, get started today."

- Speed
- Durability
- Real Gradation

Speed Amber Mill Direct provides you with the following options: All you need is 1 hour! From the time the patient sits to the time that you deliver your restoration, Amber Mill Direct speeds up your workflow. Polish Scan & Design Polish or Cementation After milling, just polish Stain / Glaze (about 9 min 10 sec) and deliver the resto-* Single crown case ration directly to the patient. Achieve excellent aesthetic results with our gradated translucency PRE-CRYSTALLIZED without any firing. Amber Mill Direct is a Lithium Disilicate-based millable glass ceramic block that requires no-crystallization, therefore, Stain / Glaze no oven. Given the shortened If your restoration fabrication time, one-hour requires more restorations are possible. characterization, simply stain / glaze it to achieve better aesthetic results. Cutting force (N) depth of cut: 0.05mm Amber[®] Mill Direct 70 Product T 60 Modify Opacity Cutting force (N) Simply bake at over 840 °C to modify the value and opacity of the restorations from 30 HT to LT. 20 10

10K

Cutting speed (rpm)

20K

30K

40K

Source: HASS R&D, Korea

Real gradation

Gradated translucency

Amber Mill Direct achieves natural translucency by applying a gradated microstructure, from the cervical to incisal/occlusal regions, without additional characterization.



shows similar translucency in cervical and incisal part to natural teeth.


*Contrast ration to natural teeth

· Enamel : 0.3~0.8 / 0.55~0.90

· Dentin: 0.6~0.95



Contrast Ratio (CR)



$$CR = \frac{Y_b}{Y_w}$$

 $CR = rac{Y_{b}}{Y_{w}}$ Yb and Yw is spectrum reflection ratio measured in black and white background. In CR, 0 means transparence and 1 means opaque.


Durability

Microstructure

Amber Mill Direct produces restorations with different microstructures that generate different strengths in the cervical and incisal regions, thus, reducing wear of the antagonist teeth.

Cervical







Edge stability

Achieve excellent marginal fit and cervical contour.





AmberMill Direct




Competitive product

Adley Lithium Odlice

Fracture strength

Amber Mill Direct produces restorations with different microstructures that generate different strengths in the cervical and incisal regions, thus, reducing wear of the antagonist teeth.





Chewing simulator

*1,000,000 cycles / 1.5 Hz / 10kg force (in pH 7.2 Water) and thermal cycling at 5-55 °C for 30s each


Test result from Chewing simulator proves superior wearout resistance in occlusal region.

Fracture strength before/after Chewing simulator





Fracture Load




Workflow

1. Scan



2. Design



3. Nest



4. Mill





5. Sprue Removal





6. Polish or Stain / Glaze (optional)





7. Cement



#11, 12, 21, 22 veneers Source : Dr Ana Petri / Oral Design Clinic

Product Q&A



- As a functional gradient block, Amber Mill Direct has different trans and strength for each area; how can we distinguish the incisal/cervical area?
- A The section where our product logo is marked on the block is the incisal area, which is more transparent, and the opposite side is the cervical area, which is more opaque. Take these points into consideration when you design your case.
- Q How is the gradated effect of your block different from other existing lithium disilicate-based glass ceramics?
- A Amber Mill Direct is uniquely designed to achieve the most natural gradation to resemble how a natural tooth gradates. We coined this unique feature as our GLD technology Gradient lithium-disilicate technology.
- Q Why does the Amber Mill Direct have a curved shape in the notch part of the holder?
- A The curved shape allows the targeted area to be reached faster allowing for low bur consumption and faster milling.
- Q Amber Mill Direct provides the option to change translucencies from HT to LT by co-firing. What is the heat treatment schedule to achieve LT?

Α	Stand-by temperature B	Closing time S	Heating rate t ₁	Firing temperature T ₁	Holding Time H ₁	Vacuum 1 V ₁₁ /V ₁₂	Vacuum 2 V ₂₁ /V ₂₂	Long-term cooling L	Cooling time
	400°C	3:00 min.	45°C	840°C	1:00 min.	450°C	840°C	690°C	-

- Q What are the pretreatment conditions used for cementation?
- A A silane for glass ceramics is applied after etching the case's inner surface for 20 seconds using 5% HF. After that, you can bond it using conventional self-adhesive resin cement.

Amber® Mill Direct





INDICATIONS





Onlays



Veneers



Anterior Crowns



Posterior Crowns * Occlusal wall thickness ≥ 2.0 mm

PRODUCT LINE-UP

Size C14 / HT Dimensions (mm)

14 × 12 × 18


pcs / Pack 5 blocks

AVAILABLE SHADES









Printed in the United States

HASS Bio America, Inc.