


A 836

Set of injector nozzles

For optimum reprocessing of lumened instruments used in dental applications.



- For use in A 105/1 upper basket or the A 315/1 module
- Pack of 22 dental injector nozzles (A 833)

EAN: 4002516204091 / Material number: 11054660 / Old Material Number: 69783600D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0.20
Gross weight in kg	0.20